

For Veterinary Use Only

*Best practices
fecal testing is a lot
easier, cheaper
and more practical
than you think!*

**StatSpin[®]
OvaTube[®]**

Now you can perform “gold standard” centrifugal flotation in your clinic and meet best practice standards for fecal parasite testing. The StatSpin OvaTube is clean and simple to use. It requires a minimum of space and can process a sample to a ready-to-read slide in record time so you can diagnose and treat in the same visit.


Best of all, OvaTube means added revenue. With as few as 10 tests/week, an OvaTube system makes strong economic sense ... and like many clinics that have started using OvaTube, you may see your fecal test volume increase to become one of your most profitable operations!

- Meets CAPC recommendations for intestinal parasite testing – accurate results with confidence
- Standardized method for all labs, large and small
- Diagnose and start treatment while patient is in-clinic – avoid lost medication revenue
- Realize substantial profits while improving your service to pet owners

StatSpin[®]



An OvaTube workcell takes only about 2 ft. of bench space, yet it can generate significant income and enhance your parasite detection program.



The new OvaTube Twist Rack simplifies the transfer of coverslips to microscope slides.



Ordering Information



Centrifugal accuracy without the complexity and mess!

Studies show that OvaTube is substantially equivalent to the standard centrifuge method, but there's less mess, less waiting and fewer transfer steps.

See the www.statspinvet.com website for the studies or call us for printed copies.

Centrifuge Compatibility

Ovatube works in most centrifuges that take 15 mL tubes ... best in low cost fixed angle models.


See the www.statspinvet.com website for centrifuge compatibility guide or call us for advice.



OvaTube Kit

SSOT-50

Includes Mixers, Tubes and Filters in a convenient dispensing box. 50 tests.



OvaTube Twist Rack

Makes twist for 3-minute settle easy, facilitates transfer of coverslips.

SSOT-10R

10-place rack with suction feet



NEW

OvaTube Take Home Collector

Lets pet owner hygienically collect sample for drop-off or next visit. Easy to use, simplifies processing at the lab.

SSOT-HC

50 collectors and 50 baggies

- ❑ **STOP** losing revenue and profits by sending samples out!
- ❑ **STOP** using a suboptimal method for zoonotic parasite detection!
- ❑ **START** using OvaTube, the safe, sensible and profitable alternative!

Contact StatSpin at 1-800-STATSPIN
(1-800-782-8774) for OvaTube pricing and terms.

www.statspinvet.com

A Division of IRIS International, Inc.

StatSpin®

60 Glacier Drive
Westwood, MA 02090-1825
Tel) 800.782.8774
Tel) 781.551.0100
Fax) 781.551.0036
www.proiris.com

64-006455-001 Rev. A